



The **ACIST CVi**® Contrast Delivery System

Simplifying **Contrast Injection**



In pursuit of the perfect image.

Simplifying Procedures

The **ACIST CVi Contrast Delivery System** – the sophisticated system that simplifies contrast injection for **ALL** your interventional and diagnostic cardiology procedures, from **small injections** in the coronary arteries, to **large volumes** in the ventricles and peripheral vasculature.






Simplifying Control

The ACIST CVi System's **innovative design** helps make it simple to operate, allowing you to generate **quality images** and to focus on what's important – your patient and the procedure.

Simplifying Patient Care

We are passionate about patient care. The ACIST CVi System has been shown to **reduce procedure time and the volume of contrast delivered** to the patient by providing precise contrast delivery.



- ① The **AngioTouch[®] hand controller** allows real-time, variable-flow control of the contrast injection rate for precise and consistent contrast administration, and has been **shown to reduce per-patient contrast dosage by up to 20%¹**.
- ② The **touch screen monitor** provides intuitive on-screen prompts for setup, adjustable injection volume and flow rate limits, contrast tracking information, and real-time readout, for continuous system and procedure monitoring.
- ③ The **built-in air column detection sensor** alerts the clinician and **stops the injection if air is detected** in the single-use patient tubing connected to the catheter.
- ④ The five-patient, **isolated contrast reservoir** with rapid automatic refill can **reduce contrast waste** and save time between cases.
- ⑤ In-line, continuous haemodynamic monitoring provides a **real-time pressure reading**, and the automated isolation manifold provides a **sterile barrier to the contrast reservoir**.

②



The ACIST CVi boasts an array of **advanced, built-in safety** features that provide **continuous** and **automated monitoring** of all critical systems functions*, and can deliver contrast with ease even through 4Fr catheters². By reducing overall procedure time the ACIST CVi helps to reduce radiation exposure³.



①



*Sensors are designed to aid the healthcare professional but are not intended as a substitute for proper technique, monitoring, and due diligence.


1. Anne G, Gruberg L, Huber A, et al. *J Invasive Cardiol.* 2004;16(7):360–362
2. Khoukaz S, Kern MJ, Bitam SR, et al. *Catheter Cardiovasc Interv.* 2001;52:393–398
3. Brosh D, Assali A, Fuchs S, et al. *Int J Cardiovasc Intervent.* 2005;7(4):183–187

Simplifying Workflow and Efficiency


The ACIST CVi system is designed and built to **streamline procedures** and deliver faster case turnaround, while minimizing the use of contrast.

- Up to 20% reduction
in contrast dosage to the patient¹
- Up to 40% reduction
in contrast and cost^{1,2,3}
- Up to 31% reduction
in procedure and setup time⁴

Lower average contrast dose¹ PER PATIENT IN MILLILITRES




Total contrast volume³ PER PATIENT IN MILLILITRES



1. Anne G, Gruberg L, Huber A, et al. *J Invasive Cardiol.* 2004;16(7):360–362
 2. Call J, Sacrinty M, Applegate R, et al. *J Invasive Cardiol.* 2006;18(10):469–474
 3. Brosh D, Assali A, Fuchs S, et al. *Int J Cardiovasc Intervent.* 2005;7(4):183–187
 4. Lehmann C, Hotelling M. *J Invasive Cardiol.* 2005;17(2):118–121


Procedural time PER PATIENT IN MINUTES

Setup & arterial⁴



Diagnostic & PCI

Fluoroscopy³



Diagnostic & PCI

The ACIST CVi Contrast Delivery System is intended to be used for controlled infusion of radiopaque contrast media for angiographic procedures.

ACIST CVi® Product Specifications

Flow Rates	
Contrast:	User-Responsive, pre-set Variable and Fixed rates from 0.8 to 40 ml/sec, in 0.10 ml/sec increments
Saline:	Fixed rate: 1.67 ml/sec
Volume	User-Responsive, pre-set limits with variable range of 0.8 to 99.9 ml, in 0.1 ml/sec increments
Pressure Limits	User defined from 13.79 to 82.76 bar
Fill Rate	Manual or automatic refill of 3 ml/sec
Rise Time	User-defined 0 to 1 sec, in 0.1 sec increments
Programme Routine Injection Modes	Cardiac: LCA, RCA, LV/Ao, and User Defined Peripheral Vascular: Pigtail, Selective, Microcatheter, and User Defined
Monitoring Sensors	Air Column Detect* Manifold Isolation on single-use, to-catheter tubing Contrast Source Empty Contrast Reservoir Refill and Contrast Source Isolation
Imaging Interface Synchronization**	Able to synchronize with most brands of X-ray imaging equipment
Injection Delay** or X-ray Delay**	0–99.9 sec
KVO Feature***	Range of 0.1 to 10 ml/min with 20 minute timeout; maximum of 200 ml of saline dispensed
Control Panel	27 cm (10.5") Colour Touch Screen
Flexible Mounting Configurations	Table Mount with adjustable arm or stationary stem Pedestal Cart
Syringe Size	100 ml
Disposable Kit Configurations	
Contrast Reservoir (5 patient):	Contrast Reservoir 5-patient cartridge manifold connection, and contrast tubing with spike
AngioTouch® Hand Controller & Tubing:	AngioTouch® hand controller, injection line tubing, and 3-way stopcock
Automated Isolation Manifold:	Integrated system with automated isolation-manifold, low-pressure tubing and saline spike, and supplied pressure transducer cartridge; kits with no transducer also provided
Power Requirements	Factory selectable: 100 to 120 VAC, 50–60 Hz, 10 A maximum or 200 to 240 VAC, 50–60 Hz, 5 A maximum

* Sensors are designed to aid the healthcare professional but are not intended as a substitute for proper technique, monitoring, and due diligence

** Available in synchronized peripheral mode

*** Available in peripheral mode

ACIST CVi and AngioTouch are registered trademarks of ACIST Medical Systems, Inc., registered in the US. ACIST Medical Systems, Inc., reserves the right to modify the specifications and features described herein, or discontinue manufacture of the product described at any time without prior notice or obligation. Please contact your authorized ACIST representative for the most current information. © 2011 ACIST Medical Systems, Inc. All Rights Reserved. P/N: 0511.322.01

For more information from ACIST Medical Systems, Inc.

Contact us in the US:

7905 Fuller Road
Eden Prairie, Minnesota 55344
Phone: (952) 995-9300
Fax: (952) 941-4648
USA Toll-free: 1-888-667-6648

Contact us in the EU:

ACIST Europe B.V.
Renier Nafzgerstraat 114
6221 KL Maastricht
The Netherlands
Phone: +31 43 354 5030
Fax: +31 43 354 5035

Visit our website:

www.acist.com