









### **EMPOWERING** PRODUCTIVITY

The EmpowerMR® injector helps you do more



With the **intuitive**, **user-friendly** interface



Experience more operational flexibility

With the compact, rotating head

# Help increase productivity

With the dual initializeand-replace syringe



# Eliminate limits on scanner proximity

With **nonmagnetic** injector-head construction

# Achieve consistent image quality

By **eliminating RF transients** with scanners up to 7T



**EmpowerMR**®

Injector System

RF=radio frequency.

The individuals who appear are for illustrative purposes only. All persons depicted are models and not real patients.





### **EMPOWERING** SAFETY

The EmpowerMR® injector features engineering designed to protect your patients

### Avoid battery failure

With continuous power supplied by the **battery-free hydraulic drive system** 



#### Help minimize air embolism risk

With the tilt lockout feature

#### **Monitor motion**

With the color-coded, illuminating handle that alerts you to injector ram motion



## Spot possible occlusions

With **pop-up messages** of overpressure during KVO

# Identify compromised renal function

With the patient GFR calculator



 $\label{eq:GFR} \textit{GFR=glomerular filtration rate; KVO=keep veins open.}$ 

The individuals who appear are for illustrative purposes only. All persons depicted are models and not real patients.



### **EMPOWERING** DECISION-MAKING



IRiSMR® promotes effective data management, providing information essential to optimizing care

# Get critical insights

By **capturing and consolidating** essential contrast and syringe-use data



#### Manage effectively

By **generating reports** to help determine budgets, control cost, and enhance workflow

# Protect your data

By **automatically archiving** injection information



# Analyze and review easily

By **exporting data** to desktop office support programs



#### Make timely decisions

With **simplified summary screens** that provide real-time data

The individuals who appear are for illustrative purposes only.





### **EMPOWERING** SUPPORT

Give your practice the advantage of services from a world leader in imaging

## Make the personal contrast choice

with agents based on your patient's unique needs

## Save time and enhance quality control

with an effective data management tool like IRiSMR®

#### Stay up to date

with peer-to-peer medical education

#### Manage your practice

with dedicated reimbursement and account management specialists



#### **ORDERING INFORMATION**

#### **Basic System Description**

#### **EmpowerMR® Injector System**

Includes: injector head with hydraulic controller, touchscreen remote with desktop mount, electro-hydraulic interconnect, communication and power cables (includes pedestal mount for injector head)

#### FastLoad™ Syringes

2 disposable 100 mL FastLoad syringes with 96" (y) low-pressure connecting tube and bottle spikes

**IRISMR Remote Viewer** 





#### Your path to **empowerment** starts here

The individuals who appear are for illustrative purposes only. All persons depicted are models and not real patients.

#### EmpowerMR® Injector System INDICATIONS:

The EmpowerMR Injector System is indicated for vascular administration of contrast media and flushing media in conjunction with magnetic resonance imaging (MRI).

Protocols have been independently developed and are not intended as medical advice. Bracco and Bracco Injeneering S.A. shall not be responsible for any physicians' reliance on these or any other protocols.

Bracco Injeneering S.A. reserves the right at any time and without notice, to change the specifications and features described herein, or to change the production or adjust the product described.

Not all products available in all global markets.

EmpowerMR is distributed by Bracco Diagnostics Inc.; 259 Prospect Plains Road, Bldg. H; Monroe Township, NJ 08831 USA; Phone: (800) 631-5245; Fax: (609) 514-2424; Customer Service: 1-877-BRACCO 9 (1-877-272-2269); Professional Services: 1-800-257-5181 (Option 2); EmpowerMR Email: empowerinfo@diag.bracco.com; Website: imaging.bracco.com/us-en

EmpowerMR is manufactured by Bracco Injeneering S.A., Avenue de Sevelin 46  $\rm CH-1004$  Lausanne, Switzerland. http://imaging.bracco.com

#### About Bracco Injeneering S.A.

Bracco Injeneering S.A., located in Lausanne, Switzerland, is the most recently added Business Unit of Bracco Imaging. Bracco Injeneering is committed to develop the best-in-class integrated injection solutions with a strong heritage in research and innovation. It provides quality solutions for state-of-the-art radiology centers, offering proven injection technology, built on Bracco Imaging's expertise.

Thanks to this strengthened product portfolio, including the CT Exprès® 3D Contrast Media Delivery System, EmpowerCTA®+ Injector System, EmpowerMR Injector System, NEXO® Contrast Management System and NEXO [DOSE]® Multi-Modality Radiation Informatics, Bracco Imaging will be able to focus on constant innovation, not only for devices, but also for software development and data management.

CT Exprès, EmpowerCTA, EmpowerMR, IRiSMR, NEXO and NEXO [DOSE] are registered trademarks of Bracco Injeneering S.A., registered in the U.S.

FastLoad is a trademark of Bracco Injeneering S.A.

© 2017 Bracco Diagnostics Inc. All Rights Reserved. 17-020217C

